

PROSAT Sterile Delta Wipes

100% nonwoven wipes presaturated with 70% IPA & 30% deionised water

PROSAT Sterile Delta Wipes are made from a 100% nonwoven polyester fibres which are hydroentangled together using no binders or additives. The resulting fabric is extremely soft, clean and strong. The wipes are very low in particles, fibres, and extractables.

PROSAT Sterile Delta Wipes are a cost-effective polyester wipe with excellent cleanliness and wet strength, the wipes can be used for many general cleaning and disinfecting applications.

The wipes are sterilized by gamma radiation and validated sterile to AAMI Guidelines so are ideal for use in higher grade pharmaceutical cleanrooms. Presaturated wipes ensure consistent saturation of each wipe independent of operator.

Presaturated wipes can increase solvent control and accountability as well as reduce VOC emissions. The wipes are provided in convenient and easy to use peel and reseal pouches.

When used as a disinfectant, the IPA wipes are efficacious against bacteria in 1 min and yeasts in 3 mins. PROSAT Delta wipes are authorised for sale in the EU and United Kingdom under the EU and GB Biocidal Products Regulation.



FEATURES	BENEFITS
100% hydroentangled polyester fabric	• Highly sorbent with good wet strength • Low in particles and fibres • Cost effective polyester wipe
No binders or additives	• Extremely soft fabric • Leaves no residue on the cleanroom surface
Presaturated wipes	• Reduces solvent usage and VOC emissions • Ensures consistent saturation of each wipe independent of the operator
Resealable pouch	• Maintains saturation levels throughout use and makes wipe removal easy
Validated sterile via gamma irradiation to a 10⁻⁶ Sterility Assurance Level	• Suitable for use in Grade A/B cleanrooms

Description	Part No.	Size	Packaging
PROSAT Sterile Delta Wipes Presaturated with 70% IPA and 30% DI Water	PS-7030IR-BPR	230 x 230mm	20 wipes per pouch; 40 pouches per case


For more information or to request a sample,
 email infoeu@contecinc.com or phone +33 (0)2 97 43 76 98.

Use biocides safely. Always read the label and product information before use.

Copyright© 2021 Contec, Inc. All rights reserved. INT188 070921

Product Information

Material	100% Polyester (PET) 
Construction	Hydroentangled
Saturant	70% IPA (USP Grade) with 30% Deionised Water
Sterility	Gamma irradiated to a validated SAL 10 ⁻⁶
Shelf life	2 years from manufacturing date
Environment	ISO 5 - 8 Grade A/B



Technical Data

Attribute (units)	Typical Value	Test Method
Basis weight; nominal (g/m ²)	67	Contec Method
Non-volatile residue, NVR		IEST-RP-CC004.3, Sec. 7.1.2
In deionized water; (g/m ²)	0.055	
In isopropanol; (g/m ²)	0.008	
Specific ions		IEST-RP-CC004.3, Sec. 7.2.2
Sodium; (ppm)	1.9	
Chloride; (ppm)	4.7	
Particles, readily releasable		
P ≥ 0.5µm; (x10 ⁶ /m ²)	39.3	IEST-RP-CC004.2, Sec. 5.1
Fibres > 100µm; (x10 ³ /m ²)	10.4	IEST-RP-CC004.2, Sec. 5.2

VOC Content

	Per case / kg	Per pouch / kg
PS-7030IR-BPR	7.07	0.18

Efficacy Information

Test	Organism	Log Red'n	Time	Test	Organism	Log Red'n	Time
EN16615	<i>E. hirae</i>	>7.47	1 min	EN16615	<i>P.aeruginosa</i>	>7.97	1 min
EN16615	<i>S. aureus</i>	>7.68	1 min	EN16615	<i>C.albicans</i>	>6.45	3 min

Packaging

Packaging Materials						
Outer bags	Low density polyethylene (LDPE) 					
Case	Corrugated fibreboard (PAP) 					
Packaging Configuration	EA/PCH	PCH/OB1	OB1/OB2	OB2/OB3	OB3/CS	EA/CS
PS-7030IR-BPR	20	1	10	1	4	800

EA = Wipe, PCH = Pouch, OB = Outer Bag 1/2/3/ CS = Case

Recycling Key

PET	
LDPE	
PAP	

Notes

- The data shown are typical values and should not be used as product specifications.
- Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- Current and/or comparison data may be available. Please contact a Contec sales representative for details.

Test Methods

- CTM - Contec Test Method
- IEST-RP-CC004.3 Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL.